

according to Regulation (EC) No 1907/2006

**Synergy D6**

Revision date: 12.07.2012

Print date: 25.06.2014

Page 1 of 4

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Synergy D6

**1.2. Relevant identified uses of the substance or mixture and uses advised against****Use of the substance/mixture**

for dental use only

**1.3. Details of the supplier of the safety data sheet**

Company name:	Coltène/Whaledent AG
Street:	Feldwiesenstrasse 20
Place:	CH-9450 Altstätten
Telephone:	+41 (71) 75 75 300
Telefax:	+41 (71) 75 75 301
e-mail:	info.ch@coltene.com
Internet:	www.coltene.com
Responsible Department:	Regulatory Affairs msds@coltene.com

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture**

Indications of danger: Xi - Irritant

R phrases:

May cause sensitisation by skin contact.

**2.2. Label elements**

Danger symbols: Xi - Irritant



Xi - Irritant

**Hazardous components which must be listed on the label**

TEGDMA

**R phrases**

43 May cause sensitisation by skin contact.

**S phrases**

24 Avoid contact with skin.

37 Wear suitable gloves.

**SECTION 3: Composition/information on ingredients****3.2. Mixtures**

according to Regulation (EC) No 1907/2006

**Synergy D6**

Revision date: 12.07.2012

Print date: 25.06.2014

Page 2 of 4

**Hazardous components**

EC No	Chemical name	Quantity
CAS No	Classification	
Index No	GHS classification	
REACH No		
216-367-7	BISGMA	
1565-94-2	Xi - Irritant R36/37/38	5 - 10 %
203-652-6	TEGDMA	1 - 5 %
109-16-0	Xi - Irritant R36/37/38-43	
	Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT SE 3; H315 H319 H317 H335	
276-957-5	UDMA	
72869-86-4	R52-53	1 - 5 %

Full text of R and H phrases: see Section 16.

**SECTION 4: First aid measures****4.1. Description of first aid measures****After contact with skin**

After contact with skin, wash immediately with: Water and soap.

**After contact with eyes**

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

**After ingestion**

Seek medical attention if problems persist.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**Water. Foam. Extinguishing powder. Carbon dioxide (CO<sub>2</sub>). Sand.**Additional information**

Extinguishing materials should be selected according to the surrounding area.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Avoid contact with skin and eyes.

**6.2. Environmental precautions**

Do not empty into drains or the aquatic environment.

**6.3. Methods and material for containment and cleaning up**

Wipe up with absorbent material (eg. cloth, fleece). Clear contaminated area thoroughly.

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

Avoid contact with skin and eyes.

Restricted to professional users. Keep out of the reach of children.

according to Regulation (EC) No 1907/2006

**Synergy D6**

Revision date: 12.07.2012

Print date: 25.06.2014

Page 3 of 4

**Further information on handling**

After use replace the closing cap immediately.

**7.2. Conditions for safe storage, including any incompatibilities****Further information on storage conditions**

Substances sensitive to light.

Recommended storage temperature: 4 - 23 °C

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters****8.2. Exposure controls****Protective and hygiene measures**

When using do not eat, drink or smoke.

**Hand protection**

Single-use gloves.

**Eye protection**

Framed glasses.

**Skin protection**

Lab apron.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state:

Paste

Colour:

off-white

**Test method****Changes in the physical state**

Density (at 23 °C):

2.0 g/cm<sup>3</sup>

Water solubility:

insoluble

(at 23 °C)

**Solubility in other solvents**

Ethanol. (partially soluble)

**SECTION 10: Stability and reactivity****10.4. Conditions to avoid**

Light.

UV-radiation/sunlight.

**SECTION 11: Toxicological information****11.1. Information on toxicological effects****Sensitising effects**

May cause sensitization in susceptible people.

**SECTION 12: Ecological information****Further information**

According to the criteria of the EC-classification and labelling "dangerous for the environment" (93/21/EEC) the material/product is not to be classified as dangerous to the environment.

according to Regulation (EC) No 1907/2006

**Synergy D6**

Revision date: 12.07.2012

Print date: 25.06.2014

Page 4 of 4

**SECTION 13: Disposal considerations****13.1. Waste treatment methods****Advice on disposal**

Can be burnt together with household waste in compliance with official regulations in contact with approved waste disposal companies and with authorities in charge.

**SECTION 14: Transport information****Land transport (ADR/RID)****14.1. UN number:** ---**Inland waterways transport (ADN)****14.1. UN number:** ---**Marine transport (IMDG)****14.1. UN number:** ---**Air transport (ICAO)****UN/ID number:** ---**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information****SECTION 16: Other information****Full text of R phrases referred to under Sections 2 and 3**

- |          |   |
|----------|---|
| 36/37/38 | Irritating to eyes, respiratory system and skin.                |
| 43       | May cause sensitisation by skin contact.                        |
| 52       | Harmful to aquatic organisms.                                   |
| 53       | May cause long-term adverse effects in the aquatic environment. |

**Full text of H statements referred to under Sections 2 and 3**

- |      |                                      |
|------|--------------------------------------|
| H315 | Causes skin irritation.              |
| H317 | May cause an allergic skin reaction. |
| H319 | Causes serious eye irritation.       |
| H335 | May cause respiratory irritation.    |

**Further Information**

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*