

Material Safety Data Sheet

MSDS-09588-0-001

Rev.: G
Product: CIDEX® OPA Solution
Issue Date: 9/27/2010

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND THE COMPANY/UNDERTAKING

Supplier:

Advanced Sterilization Products
33 Technology Drive
Irvine, CA 92618

Customer service telephone: 1-800-755-5900
Emergency telephone number: 1-877-208-6653 24 hrs
Product name: CIDEX OPA Solution
Synonyms: None

2. COMPOSITION/INFORMATION ON INGREDIENTS

The ingredients at their given percentages in this product are not considered hazardous to your health.

Components	CAS Number	Weight %	ASP - OEL Data (TWA - 8 hr)	ASP - OEL Data (STEL):	ASP - OEL Data (Ceiling Limit Value):	
ortho-Phthalaldehyde (1,2 - benzenedicarboxaldehyde)	643-79-8	<1	19 µg/m ³	148 µg/m ³	Not Determined	

3. HAZARDS IDENTIFICATION

Emergency overview: CIDEX OPA Solution has been reported to cause anaphylactic-like reactions in bladder cancer patients undergoing repeated cystoscopy. CIDEX OPA Solution should not be used to reprocess instruments for patients that have shown previous sensitivity to this solution or any of its ingredients.

Properties affecting health: May cause sensitization by repeated skin contact.

Principle routes of exposure:

Oral: Not toxic by oral route. Not anticipated to be a significant route of exposure.
Eye contact: May cause eye irritation and redness.
Skin contact: Not toxic by dermal route. May cause skin irritation. Repeated exposure may cause skin dryness or cracking. Exposure to skin may cause temporary staining.
Ingestion: Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea. Exposure can irritate and discolor the tissues of the mouth, esophagus and other tissues of the digestive tract.

Company: ADVANCED STERILIZATION PRODUCTS
Product name: CIDEX® OPA Solution

Product code: ASPOPA
CO-1003686

Inhalation: May cause irritation, including but not limited to discharge, coughing, wheezing tightness of chest and throat, difficulty breathing and stinging sensation in nose and throat, tingling of mouth and lips, headache, lose of smell and dry mouth. Symptoms are temporary and reversible.

Hazard information:

Target organ effects: None
Reproductive effects: This material is not reported to cause reproductive effects in humans.
Mutagenic effects: Not mutagenic in the Ames assay.
Sensitization: May cause sensitization.

Carcinogenicity rating:

Components	CAS Number	ASP	NTP:	IARC:	California Proposition 65 List:
<i>ortho</i> -Phthalaldehyde (1,2 – benzenedicarboxaldehyde)	643-79-8	Not Determined	Not Determined	Not Determined	Not Listed

Signs and symptoms: None.

Medical conditions aggravated by exposure: Inhalation of vapor may cause asthma-like symptoms (chest discomfort and tightness, difficulty with breathing) as well as aggravate pre-existing asthma.

4. FIRST AID MEASURES

Eye contact: In the case of contact with eyes, rinse immediately with plenty of water for 15 minutes and seek medical attention.

Ingestion: Do not induce vomiting. Rinse mouth followed by drinking a large quantity of water.

Inhalation: Move to fresh air immediately. If experiencing difficulty breathing, seek medical attention.

Skin contact: Wash contaminated areas thoroughly with soap and water.
Remove contaminated clothing and wash before re-use.
Seek medical attention if irritation develops or persists.

Protection of first-aiders: None.

Notes to physician: Probable damage to the mucosa from oral exposure may contraindicate the use of gastric lavage.

5. FIRE-FIGHTING MEASURES

Flash point (°F):	Not applicable
Flash point (°C):	Not applicable
Autoignition temperature:	Not applicable
Flammable limits in air - lower (%):	Not applicable
NFPA rating:	Health: 1 flammability: 0 reactivity: 0

Fire fighting information:

Suitable extinguishing media:	Use any extinguishing agent which is suitable for the surrounding fire
Extinguishing media which must not be used for safety reasons:	None
Specific methods:	None
Special protective equipment for firefighters:	Wear self-contained breathing apparatus for fire fighting if necessary.
Hazardous combustion products:	None
Explosivity:	None
Explosion limits:	
lower:	None
upper:	None

6. ACCIDENTAL RELEASE MEASURES

Personal precautions:	Wear eye and skin protection while handling material for clean-up. Avoid breathing vapors and/or mists. Responders should wear chemical resistant gloves and suitable body protection. Where there is a risk of product contact with shoe, chemical resistant shoe coverings are recommended.
Environmental precautions:	Do not wash down sewers or waterways.
Methods for cleaning up:	If required, neutralize by sprinkling approximately 25 grams of glycine (CAS# 56-40-6) powder per gallon of CIDEX OPA Solution spill. Thoroughly blend the glycine into the spill using mop or other tools. Allow 5 minutes contact for neutralization. Pick up and transfer to properly labeled containers. Allow neutralization to continue for 1 hour and then dispose of in accordance with all applicable federal, state, and local regulations. Clean contaminated surface thoroughly.

7. HANDLING AND STORAGE

Handling:

Company: ADVANCED STERILIZATION PRODUCTS
Product name: CIDEX® OPA Solution

Product code: ASPOPA
CO-1003686

Technical measures/precautions:

Use in well ventilated area and use with appropriate exhaust ventilation, for example a minimum of 10 air exchanges per hour or as defined by state and local regulations.

Safe handling advice:

Wear appropriate personal protection. Avoid contact with skin, eyes and clothing. Remove contaminated clothing and launder before reuse.

Storage:**Technical measures/storage conditions:**

This product should be stored between 59°F (15°C) and 86°F (30°C)

Incompatible products:

Keep containers tightly closed.

Avoid contact with strong acids and bases. Keep from contact with oxidizing materials.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Engineering controls:	Ensure adequate ventilation.
Eye protection:	Eye protection required.
Hand protection:	Chemical resistant gloves required.
Skin and body protection	Wear suitable protective clothing.
Respiratory protection:	None required.
Other/general protection:	None required.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance:	Clear
Physical state:	Liquid
Color:	Light blue
pH:	7.2-7.8
Odor:	Slight Odor
Boiling temperature (°F):	212
Boiling temperature (°C):	100
Freezing point/range (°C):	0
Freezing point/range (°F):	32
Specific gravity:	1.0003g/cc
Evaporation rate:	similar to water
Water solubility:	Soluble

10. STABILITY AND REACTIVITY

Chemical stability:	Stable under recommended storage conditions.
Hazardous polymerization:	Hazardous polymerization does not occur.
Hazardous decomposition products:	Unknown on product.
Materials to avoid:	Strong acids and strong bases. Stay away from strong oxidizing agents.
Conditions to avoid:	Extremes of temperature and direct sunlight.

Company: ADVANCED STERILIZATION PRODUCTS
Product name: CIDEX® OPA Solution

Product code: ASPOPA
CO-1003686

11. TOXICOLOGICAL INFORMATION

Acute toxicity

LD50/oral/rat (mg/kg) = >5000 mg/kg

LD50 Dermal Rabbit (mg/kg): >2000 mg/kg

Local effects

Oral: Non-toxic
Eye irritation: May cause eye irritation.
Skin irritation: May cause skin irritation.
Inhalation: Unknown on product.

Chronic toxicity

Oral: Unknown on product.
Inhalation: Unknown on product.

Dermal: Unknown on product.

Subchronic toxicity

Oral: Oral administration of o-phthalaldehyde (active ingredient) to rats for 90 days resulted in a NOEL of 5mg/kg/day.
Dermal: Unknown on product.

Specific effects

Corrosive effects: Not corrosive.
Sensitization: May elicit an allergic reaction.
Target organ effects: None
Mutagenic effects: Not mutagenic in Ames test.

Reproductive effects: Unknown on product.
Developmental effects: Oral administration of o-phthalaldehyde (active ingredient) to pregnant rats indicated that at maternally non-toxic doses (less than 10 mg/kg/day) there was no developmental effect.

Carcinogenic effects: Unknown on product.

12. ECOLOGICAL INFORMATION

Ecotoxicity:

Ecotoxicity effects: Unknown on this product.
Aquatic toxicity effects: This product is harmful to aquatic organisms.

Mobility: Unknown on product.
Persistence / degradability: Unknown on product.
Bioaccumulation: Unknown on product.
Degradation: Unknown on product.

13. DISPOSAL CONSIDERATIONS

Waste from residues / unused products: Waste disposal must be in accordance with appropriate US, Federal, State, Local and International regulations.

Contaminated packaging: Do not re-use empty containers

Company: ADVANCED STERILIZATION PRODUCTS
Product name: CIDEX® OPA Solution

Product code: ASPOPA
CO-1003686

Methods for cleaning up:

If required, neutralize by sprinkling approximately 25 grams of glycine (CAS # 56-40-6) powder per gallon of CIDEX OPA Solution spill. Thoroughly blend the glycine into the spill using mop or other tools. Allow 5 minutes contact for neutralization. Pick up and transfer to properly labeled containers. Allow neutralization to continue for 1 hour and then dispose of in accordance with all applicable federal, state, and local regulations.
Clean contaminated surface thoroughly.

14. TRANSPORT INFORMATION**DOT:**

DOT UN-No:	Not applicable
DOT shipping name:	Not regulated
Hazard class:	Not applicable
Subsidiary risk (hazard class):	Not applicable
Packing group:	Not applicable
DOT reportable quantity (lbs):	Not applicable

IMO/IMDG:

Hazard class:	Not applicable
IMDG page:	Not applicable
IMDG-labels:	Not applicable
Packing group:	Not applicable
MFAG table No.:	Not applicable
Proper shipping name:	Not regulated
UN/Id No.:	Not applicable

ADR/RID:

Hazard class/packing group:	Not applicable
Item:	Not applicable
ADR/RID-labels:	Not applicable
UN/Id No.:	Not applicable
Proper shipping name:	Not regulated
TREM-card:	Not applicable

IATA/ICAO:

Hazard class:	Not applicable
Packing group:	Not applicable
Proper shipping name:	Not regulated
ID/UN No.:	Not applicable
IATA - label:	Not applicable
ERG #	Not applicable

TDG (Canada):

Status:	Not applicable
Packing group:	Not applicable

15. REGULATORY INFORMATION**SARA (311, 312) hazard class:**

Immediate health:	None
Delayed health:	None
Fire:	None
Sudden Release of	None
Pressure Hazard:	

Company: ADVANCED STERILIZATION PRODUCTS
Product name: CIDEX® OPA Solution

Product code: ASPOPA
CO-1003686

Reactivity: None

TSCA inventory list: Listed under TSCA: Yes

WHMIS:

WHMIS trade secret: None

WHMIS hazard class: None

Canada DSL inventory list: Listed on DSL: Yes

Notes:

1. SARA = Superfund Amendments and the Reauthorization Act.
2. CERCLA = Comprehensive Environmental Response, Compensations and Liability Act.
3. FIFRA = Federal Insecticide, Fungicide and Rodenticide Act
4. TSCA = Toxic Substance Control Act
5. WHMIS = Canadian Workplace Hazardous Materials Information System
6. This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the MSDS contains all of the information required by the CPR.

16. OTHER INFORMATION

This data sheet contains changes from the previous version in section(s):

None

Additional advice None

Literary Reference: None

MSDS format: North American Format - U.S. and Canada
This Material Safety Data Sheet was prepared in accordance with OSHA 29 CFR 1910.1200.

Disclaimer:

The information and recommendations contained herein are based upon tests believed to be reliable. However, Johnson and Johnson does not guarantee their accuracy or completeness NOR SHALL ANY OF THIS INFORMATION CONSTITUTE A WARRANTY, WHETHER EXPRESSED OR IMPLIED, AS TO THE SAFETY OF THE GOODS, THE MERCHANTABILITY OF THE GOODS, OR THE FITNESS OF THE GOODS FOR A PARTICULAR PURPOSE. Adjustment to conform to actual conditions of usage maybe required. Johnson and Johnson assumes no responsibility for results obtained or for incidental or consequential damages, including lost profits arising from the use of these data. No warranty against infringement of any patent, copyright or trademark is made or implied.

End of Safety Data Sheet

Company: ADVANCED STERILIZATION PRODUCTS
Product name: CIDEX® OPA Solution

Product code: ASPOPA
CO-1003686