

# Safety Data Sheet

<b>Product Name:</b>	<b>Butler Propy Paste</b>
1220B	Fine grit, bubblegum flavor
1221MA	Medium grit, mint flavor
1221CA	Medium grit, cherry flavor
1222MA	Coarse grit, mint flavor
1222CA	Coarse grit, cherry flavor
1223MA	Plus grit, mint flavor

## 1. IDENTIFICATION OF THE MATERIAL AND SUPPLIER

Product Name: Butler Propy Paste with Phosphated Fluoride  
 Product Use: Abrasive polishing prophylaxis paste  
 Distributor Name: Sunstar Americas, Inc.  
 Address: Chicago, IL 60630 USA  
 Guelph, ON N1K 1C7 Canada  
 Emergency Telephone: 1-800-752-7869 (Illinois Poison Center), US Only  
 1-312-906-6194 Outside of the US  
 Telephone: USA 1-888-777-3101  
 Canada 1-800-265-8353

## 2. HAZARDS IDENTIFICATION

Acute & Chronic Hazards: Irritating to eyes, respiratory system and skin.  
 Minimal amounts normally ingested are not harmful. Ingestion of large amounts may cause nausea, vomiting, headaches, dizziness and abdominal distress.  
 Risk Phrases: R36/37/38: Irritating to eyes, respiratory system and skin  
 R22: Harmful if swallowed.  
 Safety Phrases: S24: Avoid contact with skin  
 S37: Wear suitable gloves  
 S62: If swallowed, do not induce vomiting; seek medical advice immediately.

## 3. COMPOSITION / INFORMATION ON INGREDIENTS

Chemical Entity (Components)	CAS No.	Proportion
Sodium Fluoride	7681-49-4	2.72% w/w

## 4. FIRST AID MEASURES

Inhalation: N/A

# Safety Data Sheet

<b>Product Name:</b>	<b>Butler Propy Paste</b>
1220B	Fine grit, bubblegum flavor
1221MA	Medium grit, mint flavor
1221CA	Medium grit, cherry flavor
1222MA	Coarse grit, mint flavor
1222CA	Coarse grit, cherry flavor
1223MA	Plus grit, mint flavor

Ingestion: If greater than normal dose is swallowed, do not induce vomiting. Drink large quantity of water, milk or several ounces of Milk of Magnesia. Seek medical attention immediately.

Skin: Flush well with water for at least 15 minutes.

Eye: Flush well with water for at least 15 minutes.

First Aid Facilities: Seek Medical assistance for further treatment, observation and support if necessary.

## 5. FIRE FIGHTING MEASURES

Extinguishing Media: Suitable extinguishing media for the surrounding fire should be used.

Specific Hazards: None known

Precautions in connection with fire: Use of self contained breathing apparatus is recommended. Wear protective clothing to prevent contact with skin and eyes.

## 6. ACCIDENTAL RELEASE MEASURES

Emergency Procedures: Collect material into acceptable containers.

Clean Up Procedure: Mop area and clean with soapy water.

## 7. HANDLING AND STORAGE

Handling: Use original container.

Storage: Store in cool, well ventilated area.

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure Limits: Not determined

Respiratory Protection: N/A

Eye Protection: Safety glasses/goggles

Hand Protection: Protective gloves

Engineering controls: Standard room ventilation

# Safety Data Sheet

<b>Product Name:</b>	<b>Butler Propy Paste</b>
1220B	Fine grit, bubblegum flavor
1221MA	Medium grit, mint flavor
1221CA	Medium grit, cherry flavor
1222MA	Coarse grit, mint flavor
1222CA	Coarse grit, cherry flavor
1223MA	Plus grit, mint flavor

## 9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance	Firm paste with color of stated flavor
Odor	Odor of stated flavor
pH	Not determined
Vapor Pressure	Not determined
Vapor Density	Not determined
Boiling Point	Not determined
Freezing / Melting Point	Not determined
Solubility in Water	Soluble except for dental abrasives
Specific Gravity	Not determined
Flash Point	N/A
Flammability	N/A

## 10. STABILITY AND REACTIVITY

Chemical Stability:	Stable under ordinary conditions of use and storage
Conditions to Avoid:	None known
Incompatible Materials:	None known
Hazardous Decomposition	None known
Products:	
Hazardous Reactions:	None known

## 11. TOXICOLOGICAL INFORMATION

Hazardous Ingredient:	Sodium Fluoride LD <sub>50</sub> (orally) RAT: 180 mg/kg
Routes of Exposure:	Refer to Section 4 of MSDS

## 12. ECOLOGICAL INFORMATION

Ecotoxicity:	No known effects
Persistence and Degradability:	No information found

## 13. DISPOSAL CONSIDERATIONS

# Safety Data Sheet

<b>Product Name:</b>	<b>Butler Proply Paste</b>
1220B	Fine grit, bubblegum flavor
1221MA	Medium grit, mint flavor
1221CA	Medium grit, cherry flavor
1222MA	Coarse grit, mint flavor
1222CA	Coarse grit, cherry flavor
1223MA	Plus grit, mint flavor

Disposal Method and Containers      Empty containers may be disposed of in waste containers.

Special Precautions                      Disposal should be carried out in compliance with federal, state and local regulations.

## 14. TRANSPORT INFORMATION

Special Precautions for User:          Not a hazardous product

Hazchem code:                              None allocated

## 15. REGULATORY INFORMATION

Precautionary                      For professional in-office use only.

Phrases:                                  Keep out of reach of children.

Poison Schedule:                      None allocated

Note:                                      Refer to applicable international, national and local regulations and provisions.

## 16. OTHER INFORMATION

Legal Disclaimer:                      Information contained in this MSDS is based on data believed to be correct. This information does not purport to be all inclusive and shall be used only as a guide. No warranty is expressed or implied regarding the accuracy of this data or the results to be obtained from the use thereof. In no event will the manufacturer or the distributor be responsible for damages of any nature whatsoever resulting from the use of or reliance upon this information.

MSDS History                              Date of Issue:  
     Replaces Issue: 8/1/2012 from Manufacturer  
     Changes made to previous issue: Addition of Distributor name.