



# SAETY DATA SHEET

## 1. Identification of Substance and Manufacturer:

**General Purpose Cleaner Concentrate Non-Ammoniated PCN: 228**

**Use Ultrasonic Cleaner**

**Manufacturer: L&R Manufacturing company, 577 Elm Street, P.O. Box 607, Kearny NJ 07032-0607 USA**

**Publication Date 3/18/2014 Rev. F**

**Product information call 201-991-5330 [www.lrultrasonics.com](http://www.lrultrasonics.com)**

**For emergencies involving a spill, leak, fire or accident contact CHEMTREC 800-424-9300 within the United States. Or (01) 703-527-3887 (USA) for International collect calls. SUPPLIER AND EMERGENCY**

## 2. HAZARDS IDENTIFICATION



**WARNING! IRRITANT! MAY BE HARMFUL IF SWALLOWED.**

**A: WARNING STATEMENT:** Good industrial hygiene practices should be used when handling this material.

**Acute Eyes:** Direct contact causes irritation, redness and possible tearing.

**Acute Skin:** Prolonged or repeated contact causes redness, and drying of the skin.

**Acute Inhalation:** Breathing high concentrations of vapors or mists may cause irritation to the respiratory tract.

**Acute Ingestion:** Not known, but expected to cause nausea and diarrhea.

**B: POTENTIAL HEALTH EFFECTS:**

**Chronic Effects:** This product does not contain any ingredient designated by IARC, NTP, ACGIH or OSHA as probable or suspected human carcinogen.

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

Component	CAS Number	Percentage
Water	7732-18-5	50 - 70
Disodium Cocoamphodiacetate	68650-39-5	15 - 25
Propylene Glycol	57-55-6	10 - 20
Disodium Cocoamphodipropionate	68604-71-7	1 - 10
Sodium Hydroxide	1310-73-2	<1

The exact concentration of composition has been withheld as a trade secret.

## 4. FIRST AID MEASURES

**Eye Exposure:** Hold eyelids open and flush with a steady, gentle stream of water for at least 15 minutes. Seek medical attention.

**Skin Exposure:** Wash thoroughly with water. If irritation or redness develops, seek medical attention.

**Inhalation:** If respiratory irritation or distress occurs, move victim to fresh air. Seek medical attention if respiratory irritation or distress continues.

**Ingestion:** Seek immediate medical attention. DO NOT INDUCE VOMITING.

**MEDICAL CONDITIONS POSSIBLY AGGRAVATED BY EXPOSURE:** Skin contact may aggravate existing skin disease.

## 5. FIRE FIGHTING MEASURES

**FIRE HAZARD DATA:**

**Suitable Extinguishing Media:** Extinguish with dry chemical, CO<sub>2</sub> or a BC/ABC extinguisher.

**Special Fire Fighting Procedures:** None known

**Unusual Fire and Explosion Hazards:** Product will burn under fire conditions. Closed containers may explode due to build up of pressure when exposed to extreme heat.

**Hazardous Decomposition Materials:** (under fire conditions) Oxides of nitrogen and carbon.

## 6. ACCIDENTAL RELEASE MEASURES

**Cleanup and Disposal of Spill:** Mop up spill and dispose of by dilution to a sanitary sewer system as permitted by local, state and federal regulations.

**Regulatory Reporting:** Not required

## 7. HANDLING AND STORAGE

**Minimum/Maximum Storage Temperatures:** 39 to 120<sup>o</sup>F. Protect from freezing. Freezing will NOT damage the product.

Product may be thawed in warm water and then mixed to assure homogeneity.

**HANDLING:** AVOID CONTACT WITH SKIN, EYES OR CLOTHING

## 8. EXPOSURE CONTROL/PERSONAL PROTECTION

**General:** These recommendations provide general guidance for handling of this product. Because specific work environments and material handling practices vary, safety procedures should be developed for each intended application. When developing safe handling procedures, do not overlook the need to clean and dispose of the material. Waste resulting from the use of this product should be handled in accordance with Section 13: Disposal Considerations.

Exposure Guidelines: Exposure limits are recommended worker breathing limits. The following limits apply to this material:

INGREDIENTS	LIMITS	INGREDIENTS	LIMITS
Disodium Cocoamphodipropionate	None established	Disodium Cocoamphodiacetate	None established
Sodium Hydroxide	2 mg/ M3 OSHA/PEL	Propylene Glycol TWA: 10 (mg/m3) from AIHA	

**Engineering Controls:** Normal room ventilation.

**Respiratory Controls:** For reasonable uses of this material, respiratory protection should not be necessary.

**Eye/Face Protection:** Safety glasses to protect from splashing.

**Skin Protection:** Rubber or plastic gloves to avoid drying and irritation to the skin.

**Work Practice Control:** Normal hygiene in the work area should be taken when working with or handling this product.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

**Physical State:** Clear Liquid

**pH:** 11.9 to 12.5 (as concentrate)

**Water Solubility:** Soluble

**Freezing Point Range:** < 32° F

**Vapor Pressure:** Not established

**Flash Point:** >142°F (60°C)

**Viscosity:** Similar to water.

**Flammability limits (vol/vol %):** Lower: No Data Upper: No Data V.O.C. (calculated): 33.5 g/l for concentrate and 3.05 g/l used when directed.

**Method:** Tag Closed Cup

**Color:** Pale Amber

**Specific Gravity:** 1.02

**Melting Point Range:** Not available

**Boiling Point:** > 212 ° F

**Vapor Density:** Not established

**Odor:** Detergent Like Odor

**Odor Threshold:** Not available

**Evaporation Rate:** Not Applicable

**Partition Coefficient; n-octanol / water:** Not available

**Decomposition Temperature:** Not available

**Auto Ignition Temperature:** Not available

## 10. STABILITY AND REACTIVITY

**Chemical Stability:** Stable

**Conditions to be avoided:** Heat and open flame

**Materials/Chemicals to be avoided:** Strong acids, strong oxidizing agents, aluminum.

**Decomposition Type: Thermal:** Oxides of Carbon and Nitrogen

**Possibility of Hazardous Reactions: WILL NOT OCCUR**

## 11. TOXICOLOGICAL INFORMATION

**Acute Eye Irritation:** Irritating to eyes.

**Acute Skin Irritation:** No test data found for product.

**Acute Dermal Toxicity:** No test data found for product.

**Acute Respiratory Irritation:** No test data found for product.

**Acute Inhalation Toxicity:** No test data found for product.

**Acute Oral Toxicity: LD50 (rat)** > 5g/kg for Disodium Cocoamphodiacetate, **LD50 (rat)** > 5ml / kg for Disodium Cocoamphodipropionate, **LD50 (rat)** = 20 g/kg for Propylene Glycol = 20 g/kg, **LD50 (mouse)** = 40 mg/kg. for sodium hydroxide.

**Chronic Toxicity:** This product does not contain any substances that are considered by OSHA, NTP, IARC OR

ACGIH to be a probable or suspected human carcinogen. No additional test data was found for this product.

## 12. ECOLOGICAL INFORMATION

**Ecotoxicological Information:** No data found.

## 13. DISPOSAL CONSIDERATIONS

This product can be disposed in sanitary sewer system, with dilution, where permitted by local, federal and state regulations.

## 14. TRANSPORTATION INFORMATION

This product is classified as a combustible liquid for domestic shipments per 49 CFR 173.150 3 (f). Do not stack cartons more than five high.

## 15. REGULATORY INFORMATION

Inventory Issues: All components of this product are listed on the U.S. TSCA, Canadian DSL, European EINECS/ELINS chemical listings

## 16. OTHER INFORMATION

**National Fire Protection Association**

Hazard Rating, NFPA	Health	Flammability	Reactivity	Special
1	2	0	--	

### MSDS CHANGES

REV	DATE	DESCRIPTION OF CHANGE
E	4/23/13	Reviewed/Updated
F	3/18/2014	SDS Format

**Disclaimer:** The information herein is given in good faith but no warranty expressed or implied is made. This SDS has been prepared by L&R Manufacturing Company.