

according to Regulation (EC) No 1907/2006

**Cool Temp natural**

Revision date: 09.07.2013

Print date: 25.06.2014

Page 1 of 5

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Cool Temp natural

**1.2. Relevant identified uses of the substance or mixture and uses advised against****Use of the substance/mixture**

for dental use only

**1.3. Details of the supplier of the safety data sheet**

Company name:	Coltene/Whaledent AG
Street:	Feldwiesenstrasse 20
Place:	CH-9450 Altstätten
Telephone:	+41 (71) 75 75 300
Telefax:	+41 (71) 75 75 301
e-mail:	info.ch@coltene.com
Internet:	www.coltene.com
Responsible Department:	Regulatory Affairs msds@coltene.com

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture**

Indications of danger: Xi - Irritant, N - Dangerous for the environment

R phrases:

Irritating to eyes and skin.

Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

**2.2. Label elements**

Danger symbols:

Xi - Irritant; N - Dangerous for the environment



Xi - Irritant



N - Dangerous for the environment

**R phrases**

36/38

Irritating to eyes and skin.

51/53

Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

**S phrases**

26

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

37

Wear suitable gloves.

**SECTION 3: Composition/information on ingredients****3.2. Mixtures**

according to Regulation (EC) No 1907/2006

**Cool Temp natural**

Revision date: 09.07.2013

Print date: 25.06.2014

Page 2 of 5

**Hazardous components**

EC No	Chemical name	Quantity
CAS No	Classification	
Index No	GHS classification	
REACH No		
	UDMA	30 - 35 %
	Xi - Irritant R36/38	
205-016-3	diallyl phthalate	10 - 15 %
131-17-9	Xn - Harmful, N - Dangerous for the environment R22-50-53	
607-086-00-4	Acute Tox. 4, Aquatic Acute 1, Aquatic Chronic 1; H302 H400 H410	
221-950-4	TMPTMA	1 - 5 %
3290-92-4	N - Dangerous for the environment R51-53	
	Aquatic Chronic 2; H411	
01-2119542176-41		

Full text of R and H phrases: see Section 16.

**SECTION 4: First aid measures****4.1. Description of first aid measures****After contact with skin**

After contact with skin, wash immediately with: Water and soap.

**After contact with eyes**

In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

**After ingestion**

Seek medical attention if problems persist.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**Water. Foam. Extinguishing powder. Carbon dioxide (CO<sub>2</sub>). Sand.**Additional information**

Extinguishing materials should be selected according to the surrounding area.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Avoid contact with skin and eyes.

**6.2. Environmental precautions**

Do not empty into drains or the aquatic environment.

**6.3. Methods and material for containment and cleaning up**

Wipe up with absorbent material (eg. cloth, fleece). Clear contaminated area thoroughly.

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

Avoid contact with skin and eyes.

Restricted to professional users. Keep out of the reach of children.

according to Regulation (EC) No 1907/2006

**Cool Temp natural**

Revision date: 09.07.2013

Print date: 25.06.2014

Page 3 of 5

**Further information on handling**

After use replace the closing cap immediately.

**7.2. Conditions for safe storage, including any incompatibilities****Further information on storage conditions**

Substances sensitive to light.

Recommended storage temperature: 4 - 23 °C

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters****Exposure limits (EH40)**

CAS No	Substance	ppm	mg/m <sup>3</sup>	fibres/ml	Category	Origin
131-17-9	Diallyl phthalate	-	5		TWA (8 h)	WEL
		-	-		STEL (15 min)	WEL

**8.2. Exposure controls****Protective and hygiene measures**

When using do not eat, drink or smoke.

**Hand protection**

Single-use gloves.

**Eye protection**

Framed glasses.

**Skin protection**

Lab apron.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state: Paste  
 Colour: off-white

**Test method****Changes in the physical state**

Density (at 23 °C): 1.7 g/cm<sup>3</sup>  
 Water solubility: insoluble  
 (at 23 °C)

**Solubility in other solvents**

Ethanol. (partially soluble)

**SECTION 10: Stability and reactivity****10.4. Conditions to avoid**

Light.  
 UV-radiation/sunlight.

**SECTION 11: Toxicological information****11.1. Information on toxicological effects**

according to Regulation (EC) No 1907/2006

**Cool Temp natural**

Revision date: 09.07.2013

Print date: 25.06.2014

Page 4 of 5

**Acute toxicity**

CAS No	Chemical name				
	Exposure routes	Method	Dose	Species	Source
131-17-9	diallyl phthalate				
	oral	ATE	500 mg/kg		
3290-92-4	TMPTMA				
	oral	LD50 mg/kg	> 5000	Rat.	
	dermal	LD50 mg/kg	> 3000	Rabbit.	

**SECTION 12: Ecological information****12.1. Toxicity**

CAS No	Chemical name					
	Aquatic toxicity	Method	Dose	h	Species	Source
3290-92-4	TMPTMA					
	Acute fish toxicity	LC50	2 mg/l	96	Oncorhynchus mykiss	
	Acute crustacea toxicity	EC50	> 9.22 mg/l	48	Daphnia magna	

**12.3. Bioaccumulative potential****Partition coefficient n-octanol/water**

CAS No	Chemical name	Log Pow
3290-92-4	TMPTMA	3.53

**SECTION 13: Disposal considerations****13.1. Waste treatment methods****Advice on disposal**

Can be burnt together with household waste in compliance with official regulations in contact with approved waste disposal companies and with authorities in charge.

**SECTION 14: Transport information****Land transport (ADR/RID)****14.1. UN number:** ---**Inland waterways transport (ADN)****14.1. UN number:** ---**Marine transport (IMDG)****14.1. UN number:** ---**Air transport (ICAO)****UN/ID number:** ---**Other applicable information**

Not a hazardous material with respect to these transportation regulations.

**SECTION 15: Regulatory information**

according to Regulation (EC) No 1907/2006

**Cool Temp natural**

Revision date: 09.07.2013

Print date: 25.06.2014

Page 5 of 5

**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information****SECTION 16: Other information****Full text of R phrases referred to under Sections 2 and 3**

22	Harmful if swallowed.
36/38	Irritating to eyes and skin.
50	Very toxic to aquatic organisms.
51	Toxic to aquatic organisms.
51/53	Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.
53	May cause long-term adverse effects in the aquatic environment.

**Full text of H statements referred to under Sections 2 and 3**

H302	Harmful if swallowed.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
H411	Toxic to aquatic life with long lasting effects.

**Further Information**

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*