

SAFETY DATA SHEET (SDS)

In accordance with Hazard Communication Standard (HCS), 29CFR 1910.1200(g), 2012 revision

SECTION I - PRODUCT / COMPANY IDENTIFICATION

Product Trade Name: SUPERCURE Core Build-Up
Synonyms: Light-Cure Composite Restorative Paste
Recommended Use: Dental Material
Part (Item) Number: 310140, 310141, 310142, 310143, 310146, 310147
Manufacturer: CENTRIX, INC. 770 River Road, Shelton, CT 06484 USA
Emergency Telephone Number: (203) 929 - 5582
Information Telephone Number: (800) 235 - 5862

SECTION II - HAZARDS IDENTIFICATION

For potential health hazards: SEE Section 11 below.
No other significant hazards to man or the environment are reasonably foreseen.

SECTION III - COMPOSITION / INFORMATION ON INGREDIENTS

Product Contents: Light-Cure Composite Restorative Paste
Hazardous Ingredients: None listed.
Exact % content withheld due to trade secret restrictions.

SECTION IV - FIRST AID MEASURES

Skin: Wash with soap and water
Ingestion: Induce vomiting if large amount of material is swallowed
Inhalation: No harm anticipated from materials in quantities packaged.
Eyes: Flush with copious amounts of water

SECTION V - FIRE-FIGHTING MEASURES

Suitable extinguishing media: Water contact acceptable
Forbidden extinguishing media: None
Special firefighting measures: None in quantities packaged
Unusual fire and explosion hazards: None
Special protection equipment: None

SECTION VI - ACCIDENTAL RELEASE MEASURES

Personal precautions: Avoid contact with skin and eyes.
Wash hands with soap and water if contact occurs
Environmental precautions: Dispose of in accordance to all federal, state & local regulations
Reclaiming methods: Collect using mechanical means and solid absorbents

SECTION VII - HANDLING AND STORAGE

Handling: Not applicable
Storage: Store at temperatures not exceeding 75°F (24°C)
Containers: Not applicable
Incompatibilities: Not applicable
Other precautions: Avoid exposure to moisture and light

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

General measures: Not required
Limits: Not determined
Respiratory protection: Not required
Eye protection: Safety glasses or goggles
Skin protection: Not required under normal handling
Hands protection: Avoid prolonged contact skin: wash with soap and water
Other equipment:

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Tooth shade paste Explosion limits: Not Applicable
Odor: Negligible Water solubility: Virtually insoluble
Boiling point: Not Applicable Evaporation rate: Not Applicable
Flash point: Not Applicable
Other Data: Not available or none.

SECTION X - STABILITY AND REACTIVITY

Stability: Generally stable
Conditions to avoid: Avoid Exposure to elevated temperatures and light
Materials to avoid: Not Applicable
Hazardous decomposition products: Not Applicable
Hazardous polymerization products: May occur if exposed to light

SECTION XI - TOXICOLOGICAL INFORMATION

Eye contact: May cause irritation if not removed promptly
Skin contact: Possibility of skin irritation if not washed off promptly, or upon prolonged or repeated exposure
Prolonged exposure: Possibility skin sensitization or irritation upon prolonged or repeated exposure
Ingestion: Small quantities are not considered harmful
Inhalation: None

SECTION XII - ECOLOGICAL INFORMATION

No ecological damage or impact anticipated upon exposure to environment in quantities as packaged.

SECTION XIII - DISPOSAL CONSIDERATIONS

Do not dispose into sink. Dispose in accordance with all federal, state and local regulations.

SECTION XIV - TRANSPORT INFORMATION

Avoid elevated temperatures and high intensity light

SECTION XV - REGULATORY INFORMATION

U.S. FDA Dental Device: Class II product
EC Medical Device Directive: Class IIa product

SECTION XVI - OTHER INFORMATION

Product handling exclusively by dental specialists per instructions for use.

CAUTION

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

Products covered by one or more of the following patents:
US 4,963,093; US 5,122,057; and US D353673.
Other U.S and International patents pending or granted.



CENTRIX INCORPORATED

770 River Road, Shelton, CT 06484 USA
(203) 929-5582 (800) 235-5862 FAX (203) 929-6804
www.centrixdental.com



CE
0470

EC REP Dental-Contact Vertriebs OHG, Im Spieß 101, DE 26506 Norden

EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street
Sydney, NSW 2000, Australia