Beautifil II.xls

MATERIAL SAFETY DATA SHEET

1.Identification of the Product and the Company

1.1 Product Name:

Beautifil II

1.2 Company /Undertaking identification

Manufacturer's Name: SHOFU INC.

Address: 11-Kamitakamatsu-Cho Fukunine Higashiyama-ku Kyoto 605-0983 ,JAPAN

Phone Number for Information: 81-75-561-1112

Signature of Preparer:

Date Prepared:

2006.9.11.

1.3 Emergency Telephone Number

81-75-561-1112

2.Composition/Information on Ingredients

Composition

%

Bis-GMA

≒ 7.5

Triethylenglycol dimethacrylate

< 5

Aluminofluoro-borosilicate glass

⇒70

 Al_2O_3

DL-Camphorquinone

Others

Hazardous ingredients are not include.

3. Hazardous Identification

3.1 Route of Entry:

Inhalation: N/A,

Skin:

may be, Ingestion: may be,

3.2 Health Hazard(Acute and Chronic)

May cause irritation to skin and eyes.

3.3 Carcinogenicity

N/A

4.First Aid Measures

4.1 Eye contact: Flash with flowing water and seek medical advice if necessary.

4.2 Skin contact: Wash with mild soap and water.

4.3 Ingestion:

Consult physician if necessary.

5.Fire-fighting Measures

5.1 Suitable Extinguishing Media:

CO2, dry chemical, foam

5.2 Extinguishing Media Should not be used :

5.3 Special Exposure Hazards Arising:

in case of fire, do not breath fumes and gases.

(Irritation gases and fumes may emit.)

5.4 Special Protective Equipment:

Wear fire protective cloth and self-contained

breathing apparatus if necessary.

6.Accidental Release Measures

6.1 Personal Precautions:

N/A

6.2 Environmental Precautions:

Dispose according to local, state and federal regulations.

6.3 Methods for Cleaning Up;

Wipe up and discard in a suitable container.

Beautifil II.xls

MATERIAL SAFETY DATA SHEET

7.Handling And Storage

7.1 Handling:

Avoid inhalation of grindings and prolonged skin contact with uncured resin.

7.2 Storage:

Store in a cool and dark area. (Keep refrigerated when not in use)

8. Exposure Controls/Personal Protection

8.1 Respiratory Protection:

Not required(use dust mask while grinding or finishing)

8.2 Hand Protection:

Latex or rubber gloves

8.3 Eye Protection:

Safety goggles (May be required)

8.4 Skin Protection:

Not required

9. Physical And Chemical Properties

9.1 Appearance/Odour/pH

Colored putty-like paste with sweet odor

9.2 Boiling Point/Boiling Range:

N/A

9.3 Melting Point/Melting Range:

N/D

9,4 Flash point:

N/A N/D

9.5 Auto flammability:

M/D

9.6 Explosive Properties:

9.7 Oxidizing Properties;

N/A

9.8 Vapour Pressure: 9.9 Relative Density:

2.25(water=1)

9,10 Solubility:

water solubility

Insoluble

10.Stability And Reactivity

10.1 Condition to Avoid:

Visible light on liquid

10.2 Materials to Avoid:

Strong oxidizing materials

10.3 Hazardous Decomposition Products:

Will not occur

11.Toxicological Information

12.Ecological Information

13.Disposal Conditions

14.Transport Information

15.Regulatory Information

16.Other Information

Handling of this product exclusively by dental specialists according to instruction for use.