

Material Safety Data Sheet PROPHY PASTE

Page 1 of 5

Section 1 – Identification

Product Name: Gelato Propphy Paste
Manufacturer: Deepak Products LLC.
5220 N.W. 72nd Avenue Miami, FL 33166
Active Ingredient: Sodium Fluoride
Information Contacts: (305)-482-9669
Toll Free: 1-877-8-DEEPAK
Emergency Phone Numbers: US & Canada 1 (877) 8-DEEPAK

Family: Propphy Paste
Product Use: Professional Dental Polishing Paste
Product #: Various

Section 2 – Hazards Identification

EMERGENCY OVERVIEW

This information is based on findings from related or similar materials.

- Will cause eye irritation.
- May cause skin irritation.
- Prolong exposure may cause respiratory system irritation.
- Ingestion of large amounts may cause nausea, vomiting, headaches, or dizziness.
- Please read entire MSDS for additional information.



Potential Health Effects, Signs and Symptoms of Exposure:

Primary Route of Entry Inhalation, skin, and ingestion.

Eye Vapors are irritating to the eyes. Splashes may cause redness of the eye.

Skin Repeated/prolonged contact may cause mild irritation of the skin. Symptoms include redness, burning, or itching.

Ingestion Swallowing small amounts during normal handling is not likely to cause harmful effects; swallowing large amounts may be harmful. Symptoms include irritation of the throat, nausea, dizziness, and upset stomach.

Inhalation If fine small particles are inhaled, irritation of the throat with a feeling of tightness in the chest may occur.

NOTE: Refer to Section 11, Toxicological Information for Details

Section 3 – Composition/Information on Ingredients

| Chemical Identity | CAS Numbers | EINECS# | INCI Name | Exposure OSHA TWA/STEL | Limits ACGIH TWA/STEL | Carcinogen IARC/NTP/OSHA | % |
|--|--|-----------|-----------------|--|--------------------------------------|-----------------------------|-------|
| Glycerol | 56-81-5 | 200-289-5 | Glycerin | (as respirable, mist) TWA = 5 mg/ m ³ | (as mist) TWA = 5 mg/ m ³ | Not listed | 20-30 |
| Sodium Fluoride | 7681-49-4 | 231-667-8 | Sodium Fluoride | 3ppm / 6ppm | 2.5 mg / m ³ | Not Listed | 2-3 |
| N/E – None Established N/R – Not Reviewed | N/DA – No Data Available N/A – Not Applicable | | | | | | |

(items in parenthesis relate to 1999/45/EC)

Sodium Fluoride: Danger Symbol – GHS06 (T), GHS07 (Xi) Hazard Statement – H301 (R25), EUH032 (R32), H315 (R38), H319 (R36)
Precautionary Statement – P102 (S2), P260 (S23), P280 (S36/37/39), P309+314 (S45)

Glycerol: Danger Symbol – GHS07 (Xi) Hazard Statement – H315 (R38), H319 (R36) Precautionary Statement – P404 (S7), P403+ 235 (S3/9)

See Section 16 for Hazard and Precautionary Statement Key.

Section 4 – First Aid Measures

First Aid for Eye Flush with water for 15 minutes, including under eyelids. Get medical help if discomfort persists.

First Aid for Skin Wash thoroughly with soap and water. Remove contaminated clothing. Get medical help if discomfort persists.

Material Safety Data Sheet PROPY PASTE

Page 2 of 5

First Aid for Ingestion If greater than normal dose is swallowed, do not induce vomiting. Drink large quantity of water or milk. Seek medical attention.

First Aid for Inhalation Move to fresh air. If having breathing difficulty, give oxygen. Seek medical attention if discomfort persists.

Section 5 – Fire Fighting Measures

| Flash Point (°F/°C) | Flammable Limit (vol%) | Auto-ignition Temperature (vol%) |
|---------------------|------------------------|----------------------------------|
| N/A | N/A | N/A |

Method:

Extinguishing Media: Suitable extinguishing media for the surrounding fire should be used.

Fire Fighting Instructions: During a fire, irritating or highly toxic gases may be generated by thermal decomposition. If potential for exposure, wear complete personal protective equipment including self contained breathing apparatus, with full face operated in pressure demand. Fight fire from a safe distance/protected location.

Unusual Hazards: Evacuate all personnel.

Section 6 – Accidental Release Measures

Spill or Release Minor spills – Clean up immediately, avoid contact with skin and eyes. Mop area and clean with soapy water

Procedures Major spills – Clear area of personnel. Restrict access to area. Avoid contact with skin and eyes, and mop / squeegee area with soapy water. Collect contaminated water and dispose of properly.

See section 8 & section 12.

Section 7 – Handling and Storage

Handling Limit all unnecessary personal contact. Avoid breathing dust / vapors. Wear appropriate PPE.

Storage Store in a cool, well ventilated area away from heat, sparks and flame. Keep containers closed when not in use. Keep in original container provided by manufacturer.

Explosion Hazard None.

Section 8 – Exposure Controls / Personal Protection

Engineering Controls General exhaust is adequate under normal operating conditions.

Personal Protective Equipment

General To identify additional Personal Protective Equipment (PPE) requirements, it is recommended that a hazard assessment in accordance with the OSHA PPE Standard (29CFR1910.132), or European Standard EN166 be conducted before using this product.

Eye/ Face Protection Safety glasses / goggles or splash shields are required when handling. Ensure eye bath is on hand.

Skin Protection No special equipment for minor exposure (i.e. when handling small quantities), but is recommended to use impermeable clothing to prevent ANY contact with this product, such as gloves, apron, boots, or whole body suit.

Respiratory Protection In case of insufficient ventilation, wear suitable respiratory equipment.

Section 9 – Physical and Chemical Properties

| Appearance | Odor & Odor Threshold | pH | Specific Gravity | Viscosity | % Volatile |
|---------------|-----------------------|-----|-----------------------------|-----------|------------|
| Colored paste | Odor of Stated Flavor | 6-8 | (H ₂ O = 1): N/A | N/A | N/A |

| Boiling Point/ Freezing Point | Decomposition Temperature | Octanol/Water Partitioning Coefficient Log Po/w | Vapor Pressure: | Vapor Density | Evaporation Rate | Ignition | Solubility In Water (20°C) |
|----------------------------------|------------------------------|---|--------------------|------------------|---------------------|----------|----------------------------------|
| N/A | N/DA | N/DA | N/A | N/A | N/A | N/A | N/A |

| Flash Point (°F/°C) | Flammable Limit (vol%) | Auto-ignition Temperature (vol%) |
|---------------------|------------------------|----------------------------------|
| N/A | N/A | N/A |

Material Safety Data Sheet PROPY PASTE

Page 3 of 5

Section 10 – Stability and Reactivity

| | |
|--|--|
| Stability: Stable | Incompatibility (Materials to Avoid): Strong Acids |
| Hazardous Decomposition Products: Toxic Fumes (HF) | Hazardous Polymerization: Will not occur |
| Conditions to Avoid: Heat and incompatible materials | |

Section 11 – Toxicological Information

| Acute Oral Toxicity | Acute Dermal Toxicity | Acute Inhalation Toxicity | Irritation – skin | Irritation – Eye |
|---------------------|-----------------------|---------------------------|-------------------|------------------|
| N/DA | N/DA | N/DA | N/DA | N/DA |
| Sensitization | Mutagenicity | Sub-chronic Toxicity | | |
| N/DA | N/DA | N/DA | | |

Section 12 – Ecological Information

Ecotoxicological Information

| Acute Toxicity to Fish | Acute Toxicity to Invertebrates | Acute Toxicity to Algae | Bioconcentration | Toxicity to Sewage Bacteria |
|------------------------|---------------------------------|-------------------------|------------------|-----------------------------|
| N/DA | N/DA | N/DA | N/DA | N/DA |

Chemical Fate Information

| | |
|-------------------------------|--|
| Biodegradability | Biodegradable. This material is not expected to significantly bioaccumulate. |
| Chemical Oxygen Demand | N/DA |

Section 13 – Disposal Considerations

Dispose of in compliance with governmental regulation (state and federal).

Dispose of container and unused contents in accordance with federal, state and local requirements. For EU Member States, please refer to any relevant Community provisions relating to waste. In their absence, it is useful to remind the user that national or regional provisions may be in force.

Section 14 – Transport Information

| | |
|--|------------------------|
| DOT (49 CFR 172) | |
| Proper Shipping Name: | Non-regulated material |
| Identification Number: | N/A |
| Marine Pollutant: | No |
| Special Provisions: | N/A |
| Emergency Response Guidebook (ERG) #: | N/A |
| IATA (DGR): | Non-regulated material |
| Proper Shipping Name: | N/A |
| Class or Division: | N/A |
| UN or ID Number: | N/A |
| Packaging Instructions: | N/A |
| Emergency Response Guidance (ICAO)#: | |
| IMO (IMDG): | |
| Proper Shipping Name: | Non-regulated material |
| Class or Division: | N/A |
| UN or ID Number: | N/A |
| Special Provisions & Stowage/Segregation: | N/A |
| Emergency Schedule (EmS)#: | N/A |
| Other Information: | N/A |

Material Safety Data Sheet PROPY PASTE

Page 4 of 5

Section 15 – Regulatory Information

US Federal Regulations

| | |
|---|---|
| Clean Air Act: HAP/ODS | This product contains the following HAP's or ODS: • NONE |
| Clean Water Act: Priority Pollutant | This product contains the following chemicals listed under the U. S. Clean Water Act Priority Pollutant and Hazardous Substance List: • None |
| FDA: Food Packaging Status | This product has not been cleared by the FDA for use in food packaging and / or other applications as an indirect food-packaging additive. |
| Occupational Safety and Health Act | This product is considered to be hazardous under the OSHA Hazard Communication Standard. It's hazards are: • None |
| RCRA | This product contains the following chemicals considered to be hazardous waste under RCRA (40 CFR 261). • None |
| SARA Title III: Section 302 (RQ) | This product contains no chemicals regulated under Section 302 as extremely hazardous substances. |
| SARA Title III: Section 302 (TPQ) | This product contains chemicals regulated under Section 304 as extremely hazardous chemicals for emergency release notification ("CERCLA" List): • Sodium Fluoride – CAS #7681-49-4 1000 # RQ |
| SARA Title III: Section 311-312: | This product is considered to be hazardous under the OSHA Hazard Communication Standard and is regulated under Section 311-312 (40 CFR 370). It's hazards are: • Sodium Fluoride – CAS# 7681-49-4 (Corrosive) |
| SARA Title III: Section 313: | This product contains the following chemicals which are subject to the reporting requirements of Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 and 40 CFR Part 372: • None |
| TSCA Section 8(b): Inventory: TSCA Significant New Use Rule: | This product does contain chemicals listed on the TSCA inventory or otherwise complies with TSCA pre-manufacture notification requirements. • Sodium Fluoride – CAS# 7681-49-4 None of the chemicals in this material have a SNUR under TSCA. |

State Regulations

| | |
|---|--|
| CA Right-to-Know Law: California No Significant Risk Rule: | Sodium Fluoride – CAS# 7681-49-4 None |
| MA Right-to-Know Law: | Sodium Fluoride – CAS# 7681-49-4 Glycerol – CAS # 56-81-5 |
| NJ Right-to-Know Law: | Sodium Fluoride – CAS# 7681-49-4 Glycerol – CAS # 56-81-5 |
| PA Right-to-Know Law: | Sodium Fluoride – CAS# 7681-49-4 Glycerol – CAS # 56-81-5 |
| FL Right-to-Know Law: | Sodium Fluoride – CAS# 7681-49-4 |
| MN Right-to-Know Law: | Sodium Fluoride – CAS# 7681-49-4 Glycerol – CAS # 56-81-5 |


International Regulations

| | |
|---|--|
| CDSL: Canadian Inventory (on Canadian Transitional List) | Sodium Fluoride – CAS# 7681-49-4 Glycerol – CAS # 56-81-5 |
|---|--|

Material Safety Data Sheet PROPY PASTE

Page 5 of 5

Labeling according to EC directives – 1272/2008 {CLP} AND 1999/45/EC (items in parenthesis relate to 1999/45/EC)

| | |
|---|---|
| European Community: | For Gelato Propy Paste (finished product): |
|  | <ul style="list-style-type: none"> DANGER SYMBOLS: GHS07 (Xi) – Warning / Irritant HAZARD STATEMENT: H316 (N/A), Causes mild skin irritation. H319 (R36), Causes serious eye irritation H335 (R37), May cause respiratory irritation. H303 (N/A), May be harmful if swallowed. PRECAUTIONARY STATEMENT: P102 (S2), Keep out of reach of children. P261 (S24), Avoid breathing dust/fume/gas mist/vapor sprays. P280 (S36/37/39), Wear protective gloves/clothing/eye protection/face protection. P309+314 (S45), If exposed or you feel unwell, get medical advice/attention. |

Section 16 – Other Information

EU Classes and Risk / Safety Phrases for Referenced Ingredients (See Section 2):

(items in parenthesis relate to 1999/45/EC)

Danger Symbols:

GHS06 (T) – Toxic



GHS07 (Xi) – Warning / Irritants

Hazard Statement: **H301 (R25)**, Toxic if swallowed. **H315 (R38)**, Causes skin irritation. **H316 (N/A)** Causes mild skin irritation. **H319 (R36)**, Causes serious eye irritation. **H335 (R37)**, May cause respiratory irritation. **H336 (R67)**, May cause drowsiness or dizziness. **EUH032 (R32)**, Contact with acid liberates very toxic gas. **H303 (N/A)**, May be harmful if swallowed.

Precautionary Statement:

P102 (S2), Keep out of reach of children. **P260 (S23)**, Do not breathe dust/fume/gas/mist/vapors/spray. **P261 (S24)**, Avoid breathing dust/fume/gas mist/vapor sprays. **P280 (S36/37/39)**, Wear protective gloves/clothing/eye protection/face protection. **P309+314 (S45)**, If exposed or you feel unwell, get medical advice/attention. **P404 (S7)**, Store in a closed container. **P403+235 (S3/9)**, Store in a well ventilated place. Keep cool.

Hazard Rating System (Pictograms)

| | |
|---|---|
| NFPA:  | HMIS:  |
|---|---|

| | |
|-------------------|------------------|
| MSDS Prepared by: | WME |
| Revision History: | 01/15/11 Initial |

The information presented herein was obtained from sources considered to be reliable. However, this information is provided without any warranty, expressed or implied, regarding its correctness or suitability for consumers intended use and/or application. For this and other reasons, we assume no responsibility and expressly disclaim liability for loss, damage or expense arising out of any way connected with the handling, storage, use or disposal of the product. This MSDS was prepared expressly for this product. Use the materials only as directed. If the product is used as a component of another product, the information contained within the MSDS may not be applicable.

Deepak components are provided on an as is basis without warranties of any kind either expressed or implied. Deepak does not warrant the use or the results of use of the materials sold on an as is basis since they are intended for remanufacturing or repackaging. It is the sole responsibility of the user to examine and determine appropriate application and regulatory requirements associated with said Deepak components.